

A Targeted Metabolomics-Based Assay Using Human Induced Pluripotent Stem Cell-Derived Cardiomyocytes Identifies Structural and Functional Cardiotoxicity Potential




J.A. Palmer¹, A.M. Smith¹, M.R. Colwell¹, V. Gryshkova², E.L.R. Donley¹, J.-P. Valentin², R.E. Burrier¹
¹Stemina Biomarker Discovery, Madison, WI, United States; ²UCB Biopharma SPRL, Investigative Toxicology, Development Science, Braine L’Alleud, Belgium

ABSTRACT

- ◆ Implementing screening assays that identify functional and structural cardiotoxicity earlier in the drug development pipeline has the potential to improve safety and the cost and time required to bring new drugs to market.
- ◆ In this study, a metabolic biomarker-based assay, Cardio *quick*Predict (Cardio^{qp}), was developed that predicts the cardiotoxicity potential of a drug based on changes in the metabolism and viability of human induced pluripotent stem cell-derived cardiomyocytes (hiPSC-CM). Assay development and testing was conducted in two phases: (1) biomarker identification and (2) targeted assay development (Palmer et al., 2020).
- ◆ In the first phase, metabolomic data from hiPSC-CM spent media following exposure to 66 drugs was used to identify biomarkers that identified both functional and structural cardiotoxicants. Four metabolites that represent different metabolic pathways (arachidonic acid, lactic acid, 2'-deoxycytidine, and thymidine) were identified as indicators of cardiotoxicity.
- ◆ In phase two, a targeted, exposure-based biomarker assay was developed that measured these metabolites and hiPSC-CM viability across an eight-point concentration curve. Metabolite-specific predictive thresholds for identifying the cardiotoxicity potential of a drug were established and optimized for balanced accuracy or sensitivity.
- ◆ When predictive thresholds were optimized for balanced accuracy, the assay predicted the cardiotoxicity potential of 81 drugs with 86% balanced accuracy, 83% sensitivity, and 90% specificity. Alternatively, optimizing the thresholds for sensitivity yields a balanced accuracy of 85%, 90% sensitivity, and 79% specificity.
- ◆ This new hiPSC-CM-based assay provides a paradigm that can identify structural and functional cardiotoxic drugs that could be used in conjunction with other endpoints to provide a more comprehensive evaluation of a drug’s cardiotoxicity potential.

METHODS



- ◆ Human iPSC-CM (iCell Cardiomyocytes², FUJIFILM Cellular Dynamics, Inc.) were plated in 96-well plates and exposed to drug for 72 hours.
- ◆ Spent media from the last 24-hour treatment period was collected for UPLC-HRMS analysis and cell viability was assessed with the CellTiter-Fluor Cell Viability Assay (Promega).
- ◆ Samples were analyzed with rapid UPLC-HRMS method optimized for the predictive metabolites .
- ◆ Each metabolite was normalized to the reference treatment samples (0.1% DMSO). Non-linear concentration-response curves were fit for each endpoint using GraphPad Prism.
- ◆ The prediction distance for each metabolite in the Cardio *quick*Predict PM was calculated at each concentration and the final PM value was determined.
- ◆ The final PM values were fit with a non-linear concentration-response curve and the “Cardiotoxicity Potential Concentration” (cTP) was predicted based on the point where the concentration response curve crossed the “Cardiotoxicity Threshold”.

RESULTS

Cardio^{qp} Accurately Predicts Functional and Structural Cardiotoxicants

Cardio^{qp} Predictions for 81 Drugs at 10×Cmax for the Prediction Model Trained for Balanced Accuracy (BAC) and Sensitivity (SEN)

NON-CARDIOTOXICANTS	Drug	BAC-Trained: Cardio ^{qp} Model	SEN-Trained: Cardio ^{qp} Model	Drug	BAC-Trained: Cardio ^{qp} Model	SEN-Trained: Cardio ^{qp} Model	FUNCTIONAL CARDIOTOXICANTS	Drug	BAC-Trained: Cardio ^{qp} Model	SEN-Trained: Cardio ^{qp} Model	STRUCTURAL CARDIOTOXICANTS	Drug	BAC-Trained: Cardio ^{qp} Model	SEN-Trained: Cardio ^{qp} Model	FUNCTIONAL & STRUCTURAL CARDIOTOXICANTS	Drug	BAC-Trained: Cardio ^{qp} Model	SEN-Trained: Cardio ^{qp} Model
	Acetylsalicylic Acid	Non	Non	Leucine	Non	Non		Astemizole	Tox	Tox		Azidothymidine	Tox	Tox		Amiodarone	Tox	Tox
	Acyclovir	Non	Non	Loratadine	Non	Non		Bepiridil	Tox	Tox		Busulfan	Tox	Tox		Amifampridine	Tox	Tox
	Adipic Acid	Non	Non	Maltol	Non	Non		Chlorpromazine	Tox	Tox		Danorubicin	Tox	Tox		Amphotericin B	Tox	Tox
	Amoxicillin	Non	Non	Methylparaben	Non	Non		Cisapride	Tox	Tox		Dexfenfluramine	Non	Non		Amsacrine	Tox	Tox
	Ascorbic Acid	Tox	Tox	Natamycin	Non	Non		Dofetilide	Tox	Tox		Dithiazanine Iodide	Tox	Tox		Anagrelide	Non	Tox
	Aspartame	Non	Non	Phenylphenol	Non	Non		Encainide	Non	Tox		Doxorubicin	Tox	Tox		Arsenic Trioxide	Tox	Tox
	Axitinib	Non	Non	Praziquantel	Non	Non		Levomethadyl Acetate	Tox	Tox		Idarubicin	Tox	Tox		Bortezomib	Tox	Tox
	Benzoic Acid	Non	Non	Ranitidine	Non	Non		Nifedipine	Tox	Tox		Imatinib	Tox	Tox		Chloroquine	Tox	Tox
	Biotin	Non	Non	Sildenafil	Non	Tox		Ondansetron	Tox	Tox		Nandrolone Decanoate	Non	Tox		Clozapine	Tox	Tox
NON-CARDIOTOXICANTS	Cetirizine	Non	Non	Sucrose	Non	Non	FUNCTIONAL CARDIOTOXICANTS	Quinidine	Tox	Tox	STRUCTURAL CARDIOTOXICANTS	Pergolide	Non	Non	FUNCTIONAL & STRUCTURAL CARDIOTOXICANTS	Crizotinib	Tox	Tox
	Cimetidine	Non	Non	Tartaric Acid	Non	Non		Sertindole	Tox	Tox		Rofecoxib	Tox	Tox		Dasatinib	Tox	Tox
	Citric Acid	Non	Non	Thiabendazole	Non	Tox		Sotalol	Tox	Tox		Rosiglitazone	Tox	Tox		Fluorouracil	Tox	Tox
	Erlotinib	Tox	Tox	Tolbutamide	Non	Tox		Terfenadine	Tox	Tox		Tegaserod	Non	Tox		Isoproterenol	Non	Non
	Gemfibrozil	Tox	Tox	Xylitol	Non	Non		Thioridazine	Tox	Tox		Telmisartan	Non	Non		Lapatinib	Tox	Tox
	Hexylresorcinol	Non	Non					Verapamil	Tox	Tox		Trastuzumab	Tox	Tox		Mitoxantrone	Tox	Tox
												Valdecixib	Non	Non		Nilotinib	Tox	Tox
																Nortriptyline	Tox	Tox
																Paclitaxel	Tox	Tox
																Sorafenib	Tox	Tox

This study evaluated a broad range of pharmaceutical compounds , demonstrating the applicability of the Cardio^{qp} assay across multiple therapeutic classes and mechanisms of cardiotoxicity. Each drug was classified as non-cardiotoxic or cardiotoxic based on the published cardiovascular effects for each drug. Cardiotoxic drugs were further stratified into three classes, functional, structural, and general (drugs that cause both structural and functional effects). The compound set contained both cardiovascular and non-cardiovascular drugs known to cause cardiotoxicity in humans, including Na⁺, K⁺ and Ca²⁺ channel blockers, antineoplastic, antiviral, cyclooxygenase-2 (COX-2) inhibitors, receptor agonists and antagonists (adreno, androgen, angiotensin II, dopamine, histamine, muscarinic, peroxisome proliferator-activated, serotonin), and tyrosine kinase inhibitors.

Predictive Metabolites are Complimentary

Balanced Accuracy-Trained Model Performance


Metabolite	Balanced Accuracy	Sensitivity	Specificity	PPV	NPV	Functional*	Structural*	Functional & Structural*
Lactic Acid	83%	65%	100%	100%	62%	80%	50%	67%
Arachidonic Acid	74%	56%	93%	94%	54%	47%	50%	67%
Thymidine	74%	52%	97%	96%	53%	33%	50%	67%
2'-Deoxycytidine	81%	65%	97%	97%	61%	67%	50%	76%
Composite Model	86%	83%	90%	93%	74%	93%	63%	90%

Sensitivity-Trained Model Performance

Metabolite	Balanced Accuracy	Sensitivity	Specificity	PPV	NPV	Functional*	Structural*	Functional & Structural*
Lactic Acid	82%	75%	90%	93%	67%	87%	56%	81%
Arachidonic Acid	78%	69%	86%	90%	61%	53%	69%	81%
Thymidine	76%	62%	90%	91%	57%	40%	50%	86%
2'-Deoxycytidine	83%	73%	93%	95%	66%	87%	56%	76%
Composite Model	85%	90%	79%	89%	82%	100%	75%	95%

*Percent of Subclass Correctly Predicted at Cardiotoxic

Metabolite Concentration-Response Curves Vary between Types and Mechanisms of Cardiotoxicity



The Metabolite Biomarkers Identified in this Study are Biologically Relevant for Cardiotoxicity

- ◆ The metabolites have key and varying roles in modulating oxidative stress and mitochondrial function and replication and are related to known mechanisms of cardiotoxicity.
- ◆ A broad range of response types was observed for each biomarker, indicating that cardiotoxic drugs may affect different components in the pathways.
- ◆ Each metabolite responds to cardiotoxicant exposure independent of changes in cell viability.

CONCLUSIONS

- The combination of 4 metabolites predicted the cardiotoxicity potential of 81 compounds with a broad range of mechanisms and therapeutic indications with ≥85% accuracy at therapeutically relevant concentrations.
- The Cardio^{qp} assay can be used to identify both functional and structural cardiotoxicants.
- Each metabolite can identify multiple mechanisms of cardiotoxicity.
- The prediction model can be used to determine the concentration a compound shows the potential to cause cardiotoxicity.
- This method can be combined with other assays or endpoints (e.g., MEA or impedance) for a comprehensive understanding of a compound’s cardiotoxicity liability.
- Future Work: Participating in the HESI multi-site study to identify chronic treatment induced cardiotoxicity using in vitro assays

Reference: Palmer, JA, et al. Toxicol Sci. 2020; 174 (2):218-240