

Defining the Reproducibility and Applicability Domain of devTOX *quickPredict*, a Human Pluripotent Stem Cell-Based Developmental Toxicity Assay



J.A. Palmer, A.M. Smith, M.R. Colwell, E.L.R. Donley, R.E. Burrier
Stemina Biomarker Discovery, Madison, WI, United States



ABSTRACT

- Assessing the accuracy, reproducibility, and applicability domain of new approach methods (NAMs) is necessary step for establishing confidence in these methods and enabling their use in a regulatory setting.
- Over 100 chemicals have been evaluated with the devTOX *quickPredict* (devTOX^{qp}) assay, which predicts the developmental toxicity potential of a chemical based on changes in human iPS cell metabolism. The assay predicted the developmental toxicity potential across this diverse set of chemicals with 87% accuracy (88% sensitivity, 86% specificity). Within individual chemical use classes (i.e., pharmaceuticals or pesticides), assay accuracy ranged from 81% to 94%, demonstrating the broad applicability of the assay.
- To further define the assay’s applicability domain, the results were separated into different pharmacological categories and performance was assessed. The assay’s sensitivity in these pharmacological categories ranged from 50% to 100% and provides insight into the assay’s biological applicability domain.
- The reproducibility of the predictive model was evaluated using independent replicates of three chemical treatments (carbamazepine, n=45; methotrexate, n=45; thalidomide, n=17) conducted by multiple scientists with multiple iPS cell lines, freeze lots and reagents over the course of 5 years. The interpolated developmental toxicity potential (dTP) values (determined using the devTOX^{qp} predictive model) were within two standard deviations of the mean for each of the chemicals, demonstrating that the assay endpoints are reproducible over time.
- These data demonstrate the importance of understanding a NAM’s biological system, its strengths, and its limitations. Taken together, these data demonstrate the accuracy, reproducibility, and broad applicability domain of the devTOX^{qp} assay and support its use as an alternative to animal models for developmental toxicity testing.

METHODS



- Human induced pluripotent stem (iPS) cells (HYRO103 or DYRO100; ATCC) were maintained in the undifferentiated state in mTeSR1 (StemCell Technologies) on Matrigel (Corning).
- Cells were plated in 96-well plates and exposed to 8 concentrations of each test article for 48 hours. Media ± test article were replaced approximately every 24 hours.
- Spent media from the last 24-hour treatment period was collected and cell viability was assessed using the CellTiter-Fluor Cell Viability Assay (Promega).
- Media samples were analyzed with UPLC-ESI-TOF-MS to determine ornithine (ORN) and cystine (CYSS) levels.
- ORN and CYSS were normalized to spiked-in internal standards (ISTD) and the median response of the reference treatment samples (0.1% DMSO). The o/c ratio was calculated by dividing the reference-normalized value of ORN by the reference-normalized value of CYSS
- Non-linear dose-response curves for the o/c ratio, ornithine and cystine response and cell viability were fit with GraphPad Prism (GraphPad Software).
- The developmental toxicity potential (dTP, o/c ratio) and toxicity potential (TP, cell viability) concentrations were predicted from the respective dose-response curves using the developmental toxicity threshold (dTT).

RESULTS

devTOX^{qp} Accurately Predicts Developmental Toxicity Potential Across a Wide Range of Chemotypes

Balanced Accuracy	Sensitivity	Specificity	PPV	NPV
87%	88%	86%	88%	86%

Developmental Toxicants (N=66)				Non-Developmental Toxicants (N=58)			
Compound	devTOX ^{qp}	Compound	devTOX ^{qp}	Compound	devTOX ^{qp}	Compound	devTOX ^{qp}
13-cis-Retinoic Acid ^A	DT	HEPP ^A	NON	Abacavir ^B	NON	Loratadine ^A	NON
5-Fluorouracil ^A	DT	Hexaconazole ^C	DT	Acebutolol ^A	NON	Methanolic ^B	NON
9-cis-Retinoic Acid ^A	DT	Hydroxyurea ^A	DT	Acetaminophen ^A	NON	Metoclopramide ^A	NON
Abacavir ^B	DT	Ketoconazole ^A	DT	Acetylcysteine ^A	NON	MEHP ^B	NON
Acetazolamide ^A	NON	Lapatinib ^B	DT	Acycloguanosine ^A	NON	Nilotinib ^B	NON
Acitretin ^A	DT	Lenalidomide ^B	DT	all-trans-Retinoic Acid ^B	DT	Novaluron ^C	NON
all-trans-Retinoic Acid ^B	DT	Lovastatin ^B	DT	Amoxicillin ^A	NON	o-Phenylphenol ^C	NON
Aminopterin ^A	DT	Methanol ^B	DT	Ascorbic Acid ^A	NON	Oseltamivir ^B	NON
Artesunate ^{A,D}	NON	Methotrexate ^A	DT	Butylparaben ^B	NON	Penicillin G ^A	NON
Atrazine ^C	DT	Methoxyacetic Acid ^B	DT	Caffeine ^B	NON	Propylene Glycol ^B	DT
Bortezomib ^A	DT	Methylmercury ^B	DT	Camphor ^A	NON	Pyriproxyfen ^B	NON
Bosentan ^B	NON	MEHP ^B	DT	Clopyralid ^C	NON	Ramelteon ^B	NON
Busulfan ^A	DT	Myclobutanil ^C	DT	Dabigatran Etxelilate ^B	NON	Resveratrol ^C	NON
Carbamazepine ^A	DT	Nilotinib ^B	DT	Desloratadine ^B	NON	Retinol ^A	NON
Chlorophacinone ^C	DT	o,p'-DDT ^C	DT	Dibutylamine ^C	DT	Saccharin ^B	NON
Cyproconazole ^C	NON	Ochratoxin A ^C	DT	Dimethyl Phthalate ^C	NON	Sitagliptin ^A	NON
Cytarabine ^A	DT	Phenytoin ^A	DT	Dimethylamine ^C	NON	Sorbitol ^A	NON
Dabigatran Etxelilate	DT	Pomalidomide ^B	DT	Diphenhydramine ^A	NON	Sotalol ^A	NON
Dimiconazole ^C	DT	Propiconazole ^C	DT	Doxylamine ^A	NON	Sucrose ^C	NON
Dinoseb ^C	NON	Pyridaben ^C	DT	Ethylene Glycol ^B	NON	Tapentadol ^B	NON
Diquat Dibromide ^C	DT	Ramelteon ^B	DT	Fipronil ^C	DT	Tegaserod ^A	NON
Doxorubicin ^A	DT	Rotenone ^C	DT	Folic Acid ^A	NON	Tetrabromobisphenol A ^C	NON
D-Penicillamine ^A	DT	Salicylic Acid ^B	DT	Glycerol ^C	NON	Thiamine ^A	NON
Endosulfan ^C	DT	Spiroamine ^C	DT	Glycolic Acid ^B	NON	Triclosan ^C	NON
Epoxiconazole ^C	DT	Thalidomide ^A	DT	Hexazinone ^C	NON	Triethylene Glycol ^C	NON
Ethylene Glycol ^B	DT	Thiacloprid ^C	NON	Imazamox ^C	NON	Triclorazole ^C	DT
Etretinate ^A	DT	ThioTEPA ^A	DT	Imazapyr ^C	NON	Zaleplon ^B	NON
Everolimus ^A	DT	Thiram ^C	DT	Isoniazid ^A	DT	Zidovudine ^B	DT
Fingolimod ^{A,D}	DT	Topiramate ^B	DT	Levothyroxine ^A	NON	Zoxamide ^C	DT
Fluazinam ^C	DT	Trifluoromethoxy ^C	DT				
Flusilazole ^C	DT	TTNPB ^C	DT				
Genistein ^C	DT	Valproic Acid ^B	DT				
Glycolic Acid ^B	DT	Warfarin ^A	NON				


Notes: HEPP: D,L-3-hydroxy-3-ethyl-3-phenylpropionamide; MEHP: Mono(2-ethylhexyl) phthalate; A: Scored against human therapeutic C_{max}; B: Scored against Daston et al., 2014⁷ or Rat LOAEL C_{max} Concentration; C: Scored against concentration threshold; D: Prediction based on active metabolite response (Artesunate: Dihydroartemisinin; Fingolimod: Fingolimod Phosphate).

Use/Industry	# DT	# Non-DT	Balanced Accuracy	Sensitivity	Specificity
Pharmaceuticals/Vitamins	39	26	90%	87%	92%
Agrochemicals	21	11	80%	86%	73%
Industrial Solvents/Additives/Byproducts	6	10	90%	100%	80%
Personal Care Products	3	8	94%	100%	88%
Food Additives/Contaminants	3	9	94%	100%	89%

- 65 pharmaceuticals were grouped based on their pharmacological mechanism of action.
- Accuracy, sensitivity, and specificity of the devTOX^{qp} assay was determined for each class.


Pharmacological Classification	# DT	# Non-DT	Accuracy	Sensitivity	Specificity
Channel Modulators	3	1	100%	100%	100%
DNA Modifiers	2	0	100%	100%	N/A
Enzyme Modulators	7	6	85%	86%	83%
Kinase Modulator	4	0	100%	100%	N/A
Nucleoside Modulator/Central Metabolite Inhibitor	6	0	83%	83%	N/A
Receptor Modulator	4	10	75%	50%	100%
Second Messenger Modulator	0	1	100%	N/A	100%
Transcription Modulator	6	1	100%	100%	100%
Other	7	7	86%	71%	100%

Example devTOX^{qp} Results for Developmental Toxicants



- Pharmaceuticals with varying mechanisms of action are accurately predicted with the devTOX^{qp} assay, with enzyme
- Sensitivity was ≥80% for all pharmacological classes except for Receptor Modulators. In many cases, the ability of the assay to accurately predict a developmentally toxic chemical is related to whether or not iPS cells express target (receptor, channel, enzyme) of interest. For example:
 - Fingolimod Phosphate (active metabolite of Fingolimod), a sphingosine 1-phosphate receptor inhibitor, decreases the o/c ratio at biologically relevant concentrations. Human iPS cell gene expression data indicate that the cells express S1PR3.
 - In contrast, the concentration where an effect is observed for Bosentan, an endothelin receptor antagonist, is >10x greater than the *in vivo* exposures. Human iPS cell gene expression data indicate low/variable expression of EDNRA and EDNRB.

devTOX^{qp} is Reproducible Over Time



Average dTP Concentrations			
	Carbamazepine	Methotrexate	Thalidomide
dTP Mean (±SD), [μM]	1.39 (±0.7)	0.06 (±0.02)	0.14 (±0.07)
dTP Range [μM]	0.31-2.77	0.03-0.12	0.05-0.26

dTP results from 45 separate experiments between 2014 and 2021 showed good reproducibility across multiple reagent lots and scientists. Over 99% (106/107) of the dTP values were within two standard deviations of the mean. On the graph, points are the individual dTP concentrations from each experiment, with a line at the mean with the standard deviation.

CONCLUSIONS

- devTOX^{qp} predicted the developmental toxicity of 115 different chemicals with 87% balanced accuracy (88% sensitivity, 86% specificity).
- devTOX^{qp} identifies the developmental toxicity potential for a wide range of mechanisms, but may not be applicable in all cases (i.e., targets not expressed in human iPS cells). Note: this is based on classification by pharmacological mechanism. The mechanism of action for developmental toxicity may be different.
- devTOX^{qp} results were reproducible across multiple reagent lots and technicians.
- Future Studies: Test all drugs included in the ICH S5(R3) Reference Compound List and Incorporate *in vitro* to *in vivo* extrapolation (IVIVE).

References:

- [1] Daston GP, et al. *Birth Defects Res B Dev Reprod Toxicol.* 2014; 101(6):423-8. doi: 10.1002/bdrb.21132; [2] Guenther MG, et al. *Cell Stem Cell.* 2010; 7(2):249-257. doi: 10.1016/j.stem.2010.06.015; [3] Exposure data from Daston et al., 2014, ICH S5(R3), or FDA Pharmacology Review for that drug. References available upon request.